
**OPERATING INSTRUCTIONS**  
FOR  
**COLD LIGHT SOURCE LED-F1/LED-F2**

power supply can change depending on availability  
 (can be either desktop OR wall plug power supply)

- 1 ... LED green (POWER)
- 2 ... LED red (TEMP.)
- 3 ... Brightness control
- 4 ... Lightguide connection
- 5 ... ESD connection
- 6 ... 24DC input
- 7 ... Type plate
- 8 ... Socket for mains cable
- 9 ... 24VDC plug connector
- 10 .On/Off switch
- 11 .Remote (optional)

# OPERATING INSTRUCTIONS FOR LED-F1/F2 COLD LIGHT SOURCE

Please read all the operating instructions completely and carefully before operating this device!

## TECHNICAL DATA

### LED-F1/F2

Model	LED-F1(R)	LED-F2(R)
Lamps	High-power LED	
Range control	min. 30:1	
Color temperature typ.	5,500 K / daylight	
Colour reproduction index	80	
Power input typ.	22W	
Power supply	24V DC +/- 5%, 1000mA Socket f. hollow pin 5.5 x 2.1 mm	
On/Off switch	✓	
Brightness control	Infinitely variable using rotary potentiometer, with Index marks 1-6	
Indicator lamps	LED RED (Temperature status) LED GREEN (Functional status)	
ESD socket	for 4mm banana plug for potential equalisation	---
Cooling	without fan by convection	with DC fan
Housing	vents on top and side	top and side closed
Remote	only with "R" models	
Lightguide connection	Photonic standard dia 15mm (galvanically isolated)	
Operating environment for normal operation	10-30°C and 30-70 % relative humidity	
Storage/transportation	0-40°C and 30-70 % relative humidity	
Dimension height x width x depth	96 x 110 x 175 mm without controls	
Weight	2500g	
Type of protection (Protection class)	III	
Conformity	CE (EMC Directive 89/336/EEC), RoHS	

(\*) LEDs are suitable for use at temperatures up to 40 °C. However, increasing temperature and humidity results in a reduction in service life. For optimal use, the ambient temperature should be 25 °C (or lower) and the relative humidity 50 % (or lower).

## Description

This new type of LED illumination device (Light Emitting Diode) has been developed for stereomicroscopy illumination purposes in industrial and laboratory settings. In purchasing this product, you have acquired a high-quality LED illumination unit for all applications requiring intensive light with minimal heat generation. LEDs are employed as the illuminant.

Some models are available as "MED" for illumination purposes in medical applications, which can be identified by the "BF" mark on the type label. When using the "MED" models, the enclosed medical power supply must be used.

Only use spareparts and accessories which are approved for use in medical applications. Contact your dealer for more information.

The following features are characteristic of this type of illumination:

- ▼ Color temperature 5,500 K (daylight!)
- ▼ Ripple and flicker free
- ▼ Long LED service life
- ▼ Quiet vibration-free operation

- ▼ Extremely compact construction

This light source has a light guide connection (4) for fibre optic products such as ring lights, goose necks, transmitted light illuminators and flexible light guides.

The modular design of the LED illumination device provides in combination with fibre-optic

- ▼ ringlights - homogeneous shadow-free illumination
- ▼ flexible lightguides – directional illumination
- ▼ transmitted light illuminators homogenous transmitted light illumination

## INDICATION FOR USE IN THE USA

If the “cTUVus” symbol is shown on the identification label, the device is built in conformance with UL2601-1. All other standards and laws required for operation for medical products (FDA, FCC...) must be met by the local dealer and demonstrated to the user. Without these proofs, the device must not be put into operation in medical purposes. In addition, note that these light sources must be equipped with a special “hospital grade” or “hospital only” power cable according to UL817/UL498 and connected to the corresponding sockets. The light source must not come into contact with the patient!!

## SYMBOLS USED ON THE DEVICE

CE Mark of conformity: Confirms conformity with EU directive.

Lights, suitable for direct attachment to normal inflammable mounting surfaces.

This symbol indicates the year of manufacture.

Optional: cTÜVus symbol: Confirms conformity with UL2601-1 if it is attached to the nameplate

This symbol on the nameplate indicates that these light sources were built for BF-type applied parts, i.e., the fiber-optics connection (flange) is isolated from the remaining device and from protective ground.

## SAFETY INFORMATION

Two warning symbols are used throughout this document:

Caution: Risk of electrical shock!

Caution: Dangerous area. Warning: Refer to accompanying documentation!

This section contains safety information which must be observed strictly when using this device. It is in your own best interest to pay attention to all warnings and safety information on the unit itself and in this manual.

### Target audience

People working with the device must read the sections which are relevant to their work. This applies especially to the chapter entitled "Safety information".

### Duty of observation with respect to product

The operator must report all operational irregularities or changes to components of relevance to safety to the responsible supervisor or the manufacturer immediately.

### **Storage location of operating instructions**

We recommend storing these operating instructions near the device to ensure quick access by operating staff when information is required.

### **Legal provisions**

National and local safety and accident prevention regulations which are in force must be strictly observed in addition to the operational guidelines issued by the operating entity.

**No technical modifications whatsoever are to be made to this device under any circumstances!**

Refer to specifications in "Technical Data" section for binding operational limits.

Do not use this unit near water or in any space with extreme humidity. In order to prevent electric shock, do not expose this appliance to rain or high humidity.

Never spill liquid on the unit or insert foreign objects into the unit. This could result in electric shock or damage to the unit.

Do not place flammable materials on or near the unit at any time. Keep the unit away from any sources of heat. The device has not been approved for operation in areas subject to explosion hazards!

The device may only be operated using the mains voltage indicated.

Never open the device or any components unless instructed expressly to do so somewhere in these instructions.

Never look directly into the LEDs or into the end of the light guide as this may result in injury or damage of your eyes.

### **Cleaning**

Before cleaning, disconnect the product from the electrical supply and let it cool off for at least 15 minutes. For cleaning, use only a cloth dampened with (soap) water or 5% ethanol or isopropyl alcohol, but never use any combustible or flammable liquids. Allow the device to dry sufficiently before startup. If fluids accidentally enter the device, unplug the power supply cable immediately and take the unit to the service location for repair. Carefully remove any dust that may have collected at the vent openings.

The use of high concentrations or an excessive cleaning can lead to scratches or discoloration of the panel!

### **Replacement parts**

Use only original spare parts. If this is not done, it can lead to personal injury and material damage. Refer to "Spare Parts" section for item numbers. If necessary contact your local dealer.

### **Responsibility**

As the ultimate legal entity, the operating institution is responsible for ensuring the proper use of the device and for providing other operators with the necessary information and it specifies the competencies required to operate the device.

Disconnect the mains cable when the unit is not being used for an extended period of time. Only use original mains cables. Route the cable so that it cannot be pinched or severed. For reasons of safety only use the grounded 3-pin mains cable.

Repairs which are not described in this document must only be carried out by authorised workshops. The manufacturer is not liable for any damage resulting from a failure to comply with the above instructions!

**Technological state**

Photonic LED illumination devices constitute state-of-the-art technology and employ recognized safety standards.

This instrument was designed and built in accordance with the following regulations and standards:

▼ EEC/89/336: Electromagnetic compatibility (CE)

**EMC warnings for "MED" models**

This chapter describes the measures required in EN60601-1-2 (electromagnetic compatibility) for ensuring electromagnetic compatibility (EMC).

Medical electrical equipment are subject to special precautionary measures with respect to EMC. Failure to observe these instructions can lead to an increased emission or a reduced noise immunity of the device or system! Portable and mobile HF communications equipment can affect medical electrical equipment!

In the medical area and application, use only accessories declared for this purpose by the manufacturer!

Do not place the device immediately next to other devices and avoid stacking with other devices. However, if the operation does require it, device and system must be closely monitored so that normal operation in this arrangement can be checked!

Use only the supplied power cables (max. 2.5 m)! Longer and different power cables can lead to an increased emission or a reduced noise immunity of the device or system!

**Table 201: Guidance and manufacturer's declaration – electromagnetic emission**

This device is intended for use in the electromagnetic environment specified below. The customer or the user of this device should assure that it is used in such an environment.		
<b>Emission test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>
RF emissions CISPR 11	Group 1	This device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.  This device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class C	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	
RF emissions CISPR 15	Complies	

**Table 202: Guidance and manufacturer's declaration – electromagnetic immunity**

This device is intended for use in the electromagnetic environment specified below. The customer or the user of this device should assure that it is used in such an environment.		
<b>Immunity test</b>	<b>IEC 60601 test level and Compliance level</b>	<b>Electromagnetic environment - guidance</b>
Electrostatic discharge IEC 61000-4-2	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines (not applicable)	Mains power quality should be that of a typical commercial or hospital environment.

Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% U <sub>T</sub> (> 95% dip in U <sub>T</sub> ) for ½ cycle 40% U <sub>T</sub> (60% dip in U <sub>T</sub> ) for 5 cycles 70% U <sub>T</sub> (30% dip in U <sub>T</sub> ) for 25 cycles < 5% U <sub>T</sub> (> 95% dip in U <sub>T</sub> ) for 5s	Mains power quality should be that of a typical commercial or hospital environment. If the use of this device requires continued operation during power mains interrupts, it is recommended that this device be powered from an uninterruptible power supply or a battery.
Power frequency (50 Hz/60 Hz) magnetic field IEC 61000-4-8	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U <sub>T</sub> is the a.c. mains voltage prior to application of the test level.		

**Table 204: Guidance and manufacturer's declaration – electromagnetic immunity**

This device is intended for use in the electromagnetic environment specified below. The customer or the user of this device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter: <b>Recommended separation distance</b>			
Conducted RF IEC 61000-4-6	3 150 kHz to 80 MHz	V <sub>eff</sub> 3 □ V1 in V	$d = \left( \frac{3,5}{V1} \right) * \sqrt{P}$
Radiated RF IEC 61000-4-3	3 80 MHz to 2.5 GHz	V/m 3 □ E1 in V/m	$d = \left( \frac{3,5}{E1} \right) * \sqrt{P}$ $d = \left( \frac{7}{E1} \right) * \sqrt{P}$ 80 MHz to 800 MHz      800 MHz to 2.5 GHz
Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (a), should be less than the compliance level in each frequency range (b). Interference may occur in the vicinity of equipment marked with the following symbol:			
At 80 MHz and 800 MHz, the higher frequency range applies. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
(a) Field strengths from fixed transmitters, such as base stations for radio (cellular /cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which this device is used exceeds the applicable RF compliance level above, this device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocation this device .			
(b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.			

**Table 206: Recommended separation distances between portable and mobile RF communications equipment and this device**

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter in Watt [W]	Separation distance according to frequency of transmitter in meters [m]		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = \left( \frac{3,5}{V1} \right) * \sqrt{P}$	$d = \left( \frac{3,5}{E1} \right) * \sqrt{P}$	$d = \left( \frac{7}{E1} \right) * \sqrt{P}$

0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters [m] can be determined using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer.

At 80 MHz and 800 MHz, the higher frequency range applies. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## INSTALLATION AND CONNECTION

Remove the individual system components from their packaging and place them onto a sufficiently large horizontal surface.

The device may only be operated at the rated mains voltage. Only connect the unit to grounded sockets.

Keep air vents free!

### **The cold light source unit must only be operated from the power mains using the power supply provided!**

Using the hollow pin plug, insert the power pack into the socket of the cold light lamp (6), providing the desktop power pack with a supply voltage. If necessary, use the ESD socket (5) for potential equalisation.


The connection from the ESD socket to the light guide (via the light guide connection) is of low impedance. In order to restrict the flow of current in the event of a malfunction, protecting the user and the equipment, it may be necessary to make provision for an additional protective resistor in the ESD connection cable.

Insert the requisite fibre optic product, e.g. ring light, gooseneck or transmitted light illuminator onto the light guide connection (4).

Never look at the lighting devices or into the end of the light guide as this may result in damage to your eyes.

### **Installation of the optional counter weight (possible only with F1(R))**

When working with goosenecks lightguides an counter weight may be necessary. Loosen the 2 screws of the cover, insert the counter weight and fix it with the 2 screws again. Store the cover on a safe place for later user.

## OPERATION AND CONTROLS

The LED cold light source can be switched on with (10).

Following a short self test, indicated by the two indicator LEDs (1 and 2) blinking, the cold light source is ready for operation. This state is indicated by the green LED (1) remaining illuminated.

### Brightness control

The brightness can be varied by turning the button (3). Please note that increasing the brightness reduces service life. Therefore, only ever work with the brightness that you really actually need.

The current brightness setting can be read on a scale. This setting is stored when the device is switched off using the button and is available when the device is switched on the next time.

### Temperature monitoring

The red LED (2) indicates if the operating temperature is too high as a result of environmental influences or a less than ideal placement of the device, resulting in a high ambient temperature. The stages are as follows:

The operating temperature approaches the permissible limit:

The red LED (2) on the front panel begins to pulsate.

If it is possible to lower the ambient temperature by taking suitable action, the LED goes off.

If this indicator is overlooked by the user, the LED starts to blink slowly and the brightness of the cold light source is reduced by approximately 30%. This enables the user to keep doing their work while the device is cooling down. If the specified operating temperature can be reached in this way, the LED control lamp (2) goes off and the LED cold light source returns to full power.

If both indicators are overlooked and/or cooling down is insufficient, the device is deactivated until the permissible operating temperature is reached. The red LED control lamp (2) is illuminated.

### Fault message

If one or both LED control lamps (1 and 2) flash rapidly during operation the device is defective. In this situation, disconnect the device from the mains, let it cool down and restart it.

If the problem still persists and one or both LEDs start blinking rapidly again, remove the device from service and contact your service department or your dealer.

Never open the device or its components.



## Service life

LED chips typically have a service life in excess of 100,000 hours. The actual service life depends to a great extent upon the ambient temperature and the operating current of the LEDs. A higher operating current produces not only more light but also more heat, which reduces the life span of the LEDs accordingly.

As a rule, the brightness of LEDs also decreases over the life cycle.

To increase the service life of the LEDs, the brightness should only ever be set to the actual level required. In this way, service life can extend to over 50,000 hours.

## Overview of LED indicator lamps

LED status		Operational state of the device	
LED state/lighting behaviour		 <b>GREEN</b> POWER	 <b>RED</b> TEMP.
Inactive		Device is disconnected from power source and therefore inactive.	Device working within permissible operating temperature.
Enabled		Device enabled	Device is overheating and switches off until it has cooled down to the specified operating temperature.
Rapid blinking		Fault message, malfunction	Fault message, malfunction
Slow blinking		Blinking during self test at device start-up. Insufficient or excessive voltage. (only use power pack supplied)	The light output of the device is reduced by approximately 30% in order for it to cool down to the permissible operating temperature (early warning level 2).
Pulsating lamp		Device is in standby mode	Operating temperature is approaching the permissible limit (early warning level).

## Spare parts and accessories

Counter weight for F1(R)	633-30-003
--------------------------	------------

Mains cable 3-pin IEC320 for:	
EU	ELN-3.001
CH	ELN-3.002
J/USA/CAN	ELN-3.003
UK	ELN-3.004
AUS	ELN-3.006
IT	ELN-3.009
USA Hospital Grade	ELN-3.010


The development of this product is undergoing continual refinement. Although we make every effort to keep these operating instructions as current as possible, we reserve the right to modify device specifications at any time without prior notification. The operating instructions are prepared subject to possible errors and misprints.

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Photonic Optische Geräte GmbH&CoKG