

F5000 MED

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Operating manual

Intended use

This product has been developed and is intended for illumination purposes in the sector microscopy and endoscopy for industry, laboratory and medical use. The applicable standards and guidelines are listed in the Chapter "Technical data". Any other use is strictly forbidden.

This product is made for use in hospitals, clinics and in doctor's offices and may be operated only by trained and qualified person. The owner is responsible for the qualification of the operator and for the compliance of all regulations and standards.

Typical applications include the provision of light for: surgery microscopes, headmounted lights and loupes, surgery endoscopes.

This user guide describes only the handling of this light source and NOT of the applicators like endoscopes or microscopes. In addition, all user guides regarding other components used (e.g. endoscope, microscope) must be observed! If multiple medical devices are combined to a system, the manufacturer of this system is responsible for the system's CE declaration.

Possible misuse


Never use this product in ophthalmology or in the vicinity of the eyes during other procedures. Misuse can lead to serious injuries to eyes.

Never use this product for intentional heating, only for illumination.

The manufacturer / distributor is responsible for and must check the compatibility of all used devices of the system. Non-observance can lead to serious injuries.

Safety Instructions

The following symbols are used in this operating manual or on the device:

	Applied part of typ CF (cardial floating)		Year of fabrication
	Manufacturer		Conformity to EC regulation
	Observe operating instruction		Protective earth
	potential equalisation		Notice: no dispose in household waste
SN	Serial number		Alternating current (AC)
	Storage and transport atmospheric pressure (min max)	REF	Article number
	Storage and transport humidity (min max)		Storage and transport temperature (min max)
	Warning of a danger area! Caution, please refer to documentation!		Warning of dangerous electrical voltage!
	Warning of hot surfaces!		ATTENTION: Dangerous optical radiation! Do not look into the light!

This chapter deals with safety instructions that must be complied with when using the device. In your own interests, please comply with all safety instructions that are listed on the device and in this manual.

The user guide must be read before using the device. It is recommended that the user guide is stored close to the device to ensure that it can be read immediately by the user if necessary. The user must inform a responsible person or the manufacturer immediately if there are any changes in the operating behaviour or to any safety-relevant device components. The operator is responsible for the correct use of the product and for the training of the user. The generally-applicable legal national safety and accident prevention regulations must be complied with in addition to the supplementary user guidelines of the operator. The manufacturer is not liable for any damages due to incorrect operation or non-compliance with the danger instructions!

Never open the device or its components, unless specifically indicated in this manual, as you could otherwise be at risk of physical injury or damage. If necessary, contact your retailer. Technical modifications to the device are not permitted under any circumstances and repairs should only be carried out by the manufacturer or authorised service points. Only original spare parts must be used. Only use original spare parts approved by the manufacturer of this product.

Do not operate the device in the vicinity of water or in rooms with extreme humidity. Do not place any liquids on or near the device. Do not spill any liquids onto the device and do not push any foreign bodies into the product, as you could be electrocuted or cause damage to the device. Do not position any flammable materials and gases (e.g. anaesthetic gas, oxygen, disinfectant,...) in the direct vicinity of the device. The device must be kept away from all sources of heat and is not approved for use in explosion hazard areas! The device may only be used in closed rooms.

Never allow the device to function unsupervised! The device may only be operated with the specified supply voltage. To reduce the risk of electric shock this device must only be used in mains power supplies with protective earth. Only use as much brightness as required for working. A higher brightness setting can irritate the eyes, stresses the object and also reduces the service life of the device.

Always keep the ventilation openings of the device free to avoid causing any overheating. Ensure that the free air circulation is not hindered by anything and that there is a sufficiently large gap between the ventilation openings and other objects. Keep the environment of the device clean (above all, dust- and lint-free), so that any device fans cannot suck in any dirt.

ATTENTION: dangerous optical radiation, do not look into the light, as it may injure your eyes

The device emits light with a high intensity! Never look directly into the light, into any light-emitting device openings or into the light guide end, as this will risk injury to the eyes. The emitted light is absorbed by materials and converted into heat. Never cover the light emission aperture of the device or the light guide outlet end with the hand, any other parts of the body or any materials. To illuminate objects correctly, always ensure there is an appropriate gap between the object and the light guide outlet end, and an appropriate brightness level to avoid thermal damage or fire. Never leave the light guide end lying around freely.

Use only fused or pressed light guides. Traditionally-bonded fibres can burn and thus destroy light guides and light source. Use only medically-approved light guides.

Note that both end pieces of the light guide, but in particular the device-side end, can become very hot. Therefore, do not touch it immediately after removal from the light source. Let the light guide cool down in the switched-off light source for a sufficient length of time before removing it.

This product, including accessories, may only be used for medical applications by a doctor or a health-care professional under the supervision of a doctor. The operator must be adequately trained. This

user guide does not contain any explanation of endoscopic or microscopic procedures – it merely describes the use of this light source. In addition, all user guides regarding other components used (e.g. endoscope, microscope) must be observed.

Never use this light source for eye surgery or surgical procedures which require direct illumination of the eye.

Note that when used with an endoscope, the latter may also emit light from its light outlet, which leads to an increase in the temperature of the outlet, of the coupling point between the light guide and endoscope and the illuminated surface. You should therefore never deploy more light than is absolutely necessary and keep the greatest possible distance between the endoscope end and the surfaces. The light-source-side light-guide connection is Type CF-insulated, allowing the safe use of energy-powered F type endoscopes.

Check before each use and after each setting change if a real live image (rather than a permanently stored one) is present, with the correct image orientation. Before each use, it must also be ensured that all parts and accessories meet the criteria given in the instruction for use of all devices, so as to ensure safe use.

Select the active light-guide diameter so that it is at most as large as the inlet diameter of the endoscope. Ideally, it will be a little smaller. If it is larger than the endoscope input, the excess light is absorbed by the connecting parts, which leads to very high temperatures at the coupling point and can subsequently cause severe burns and coagulation. Be sure to observe the user guide for the light guide, the endoscope and all other components used.

Likewise, the distal endoscope end-piece may become very warm and may also lead to the above clinical effects (tissue damage, coagulation, burns).

Before using a defibrillator this device must be removed from the surgical area.

Before the theoretical end of a bulb's lifespan, you should reckon on bulb failure, and therefore have a second functioning light source at the ready in the case of urgent and life-threatening operations. Bear in mind that in the event of a power failure, no light will be available and, as required, ensure that the light source has an uninterruptible power supply (UPS). Bulb replacement is done by the manufacturer alone, by contacting the local distributor.

When connected to an endoscope or microscope, light can emerge with high-radiation energy, which can lead to an increase in temperature of the irradiated area (tissue) and of the light outlet end-piece. Read and observe the safety and operating instructions provided by the microscope or endoscope manufacturer! If the light source is used together with other medical electrical equipment or any other part of an endoscopic device, interactions and potential hazards may arise!

If this product is combined with other medical products (e.g. endoscope, microscope, headlight...), the result is a medical system that has to be tested and approved by the user or distributor as a whole. Ensure that all components are compatible with each other and check this as required with the distributor.

The light source and light guide must not come into contact with the patient. Place the light guide so that no one is put at risk (fall, strangulation, ...) and so that no objects may be caused to fall.

EMC Warnings

This chapter describes the measures required in EN60601-1-2 (electromagnetic compatibility) for ensuring electromagnetic compatibility (EMC).

Medical electrical equipment are subject to special precautionary measures with respect to EMC. Failure to observe these instructions can lead to an increased emission or a reduced noise immunity of the device or system! Portable and mobile HF communications equipment can affect medical electrical equipment!

Use only accessories declared by the manufacturer as being for the given medical area and

application.

Do not place the device immediately next to other devices and avoid stacking with other devices. However, if the operation does indeed require this, the device and system must be closely monitored so that normal function can be assured whilst in this arrangement.

Use only the supplied power cables (max. 2.5m). Longer or different power cables may lead to increased emissions or reduced noise immunity of the device or system.

Table 201: Guidance and manufacturer's declaration – electromagnetic emission

This device is intended for use in the electromagnetic environment specified below. The customer or the user of this device should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	This device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	This device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	
RF emissions CISPR 15	Complies	

Table 202: Guidance and manufacturer's declaration – electromagnetic immunity

This device is intended for use in the electromagnetic environment specified below. The customer or the user of this device should assure that it is used in such an environment.		
Immunity test	IEC 60601 test level and Compliance level	Electromagnetic environment - guidance
Electrostatic discharge IEC 61000-4-2	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines (not applicable)	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% U_T (> 95% dip in U_T) for ½ cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles < 5% U_T (> 95% dip in U_T) for 5s	Mains power quality should be that of a typical commercial or hospital environment. If the use of this device requires continued operation during power mains interrupts, it is recommended that this device be powered from an uninterruptible power supply or a battery.
Power frequency (50 Hz/60 Hz) magnetic field IEC 61000-4-8	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to application of the test level.		

Table 204: Guidance and manufacturer's declaration – electromagnetic immunity

This device is intended for use in the electromagnetic environment specified below. The customer or the user of this device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter: Recommended separation distance			

Conducted RF IEC 61000-4-6	3 V _{eff} 150 kHz to 80 MHz	3 □ V1 in V	$d = \left(\frac{3,5}{V1}\right) * \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 □ E1 in V/m	$d = \left(\frac{3,5}{E1}\right) * \sqrt{P}$ $d = \left(\frac{7}{E1}\right) * \sqrt{P}$ 80 MHz to 800 MHz 800 MHz to 2.5 GHz

Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (b). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (a), should be less than the compliance level in each frequency range (b).
Interference may occur in the vicinity of equipment marked with the following symbol:
At 80 MHz and 800 MHz, the higher frequency range applies. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

(a) Field strengths from fixed transmitters, such as base stations for radio (cellular /cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which this device is used exceeds the applicable RF compliance level above, this device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocation this device .
(b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

Table 206: Recommended separation distances between portable and mobile RF communications equipment and this device

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter in Watt [W]	Separation distance according to frequency of transmitter in meters [m]		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = \left(\frac{3,5}{V1}\right) * \sqrt{P}$	$d = \left(\frac{3,5}{E1}\right) * \sqrt{P}$	$d = \left(\frac{7}{E1}\right) * \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters [m] can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer.
At 80 MHz and 800 MHz, the higher frequency range applies. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Operation

The safety instructions must be observed.

Unpack the light source and remove any transport protection. Place the light source on a flat, stable, vibration-free surface, ensuring that it is a sufficient distance from other objects. Save the packaging used for subsequent storage and shipping.

Check whether the device-side light-guide end-piece is clean, as impurities can lead to burning,

and subsequently to the destruction of the light source and the light guide. Due to the high luminous intensity, bonded light guides may be destroyed. Ensure that hot-fused or pressed light guides are used.

Before each use:

- Check that the light source, the light guide, the power cord and all the other parts to be used are clean, functional and free of defects, in particular the electrical and mechanical connections. Only use power cords with protective earthing and operate the unit only on mains power supplies with an intact protective earth.
- Check that the correct adapter is used for the light guide.
- Check the secure fit of the adapter and the light guide in the light source.
- Check that the potential compensation connector is plugged in (10), insofar as this is required by the operator.
- Connect the light source to the power grid only after all other components have been installed and tested.
- Check whether the light source still has sufficient brightness for the planned activities.

Now insert the light-source-side end-piece of the light guide into the designated light-source connection (6) and the other end into the endoscope or microscope.

If the system is now fully installed and tested, you can now energise the light source (8) and switch it on.

Never look into the light beam and maintain a safe distance between the light guide outlet and other objects. All safety instructions must be followed in this respect. Check the temperatures of all parts of the system (e.g. endoscope, light guide, ...), especially those parts that may come into contact with patients, in order to avoid injury.

Setting-up/removal of light guide:

Note that both end pieces of the light guide, but in particular the device-side end, can become very hot. Therefore, do not touch it immediately after removal from the light source. Let the light guide cool down in the switched-off light source for a sufficient length of time before removing it.

If the correct adapter (6) is installed, the light guide need only gently be pushed into the light-source adapter until it stops. The light guide is automatically clamped. To remove it, the adapter must be slightly pressed, thereby releasing the clamp and allowing the light guide to be removed.

If the light source is turned on but no light guide connected, the bulbs remain dark to avoid glare.

Switching on:

First, switch on the light source using the main switch (7) at the rear. The LED above the standby button (1) will now display a faint glow, signifying standby mode. You can now set the light source into operation by briefly pressing the standby button on the front (1).

Switching off:

The light source can be switched completely off using the main switch (7) at the rear, or may be switched to standby mode by pressing and holding the button (1). It is recommended that the unit is always switched off in order to save energy and to increase the lifespan of the device.

Disconnection from power supply:

The unit is completely switched off by depressing the main switch (7). Disconnection from the power supply is done by pulling the power plug from the socket (8) at the rear. The device should therefore always be positioned so that there is easy access to this power-disconnection equipment.

Brightness control:

The brightness can be controlled by operating the touch wheel (3). To do so, run a finger along the annular groove clockwise to increase the brightness (or anti-clockwise to reduce brightness). Flashing LEDs along the touch wheel show the approximate current brightness. A more accurate indication is given by the display (5), which indicates the current brightness in 100 degrees.

The brightness can also be gradually increased (+) or decreased (-) using (4).

Similarly, the preset buttons (2) can be used for brightness control. If you depress a key, the stored

brightness is adjusted and a LED indicator above the key shows the activation of the preset. Any values within the permissible range can be stored. To do this, set the desired value by means of the touch wheel (3) or + | - keys (4). Then hold down the desired preset button (2) until the LED above it blinks briefly. This value is now stored under the selected preset key.

Display:

The OLED display (5) indicates the current brightness value in 100 degrees, with 100 being the maximum brightness and 1 the minimum. 0 means that the LED is off, but that the power is still on. If the light flashes, an overheating fault or fan failure has occurred (see Troubleshooting).

Light-guide adapter:


Change the adapter only when the light source has been switched off and detached from the power source. Note that the adapter and light guide may be hot after operation; therefore, let the parts cool down sufficiently before changing.

Each type of light guide requires its own adapter. For this reason, the adapter may be easily exchanged for another.

Remove the light guide and leave all parts to cool off sufficiently. Turn the adapter (6) anticlockwise slightly and then pull it out.

Now insert the new adapter (6) into the hole and twist it until it clicks clockwise until it clicks.

Options (11), (12) and (13):

Are only for service purposes and are covered at the factory. Do not remove the covers!

Troubleshooting

This chapter helps you to troubleshoot. If the error cannot be found, the device must be sent in to the service team. The device may in no circumstances be operated if a fault exists or is suspected.

Fault	Remedy
No light	<ul style="list-style-type: none"> • Is a power cord plugged in which has been plugged into a working power outlet? • Is the rear power switch (7) turned on? • Has the front button (1) been switched and is the associated LED indicator light on brightly and is the display on? • Is there a light guide connected? • Are the fuses (9) in working order? • Is the display flashing? If so, then there is a temperature or fan problem.
Display flashing	<p>Temperature or fan problem</p> <ul style="list-style-type: none"> • Turn off the light source (7) and let it cool, then it can be switched

	<p>on again. Ensure sufficient cooling of the device.</p> <ul style="list-style-type: none"> • If the display is still flashing following the above, it is likely that there is a fan problem. In this case the device must be sent in to the service team.
Brightness decreases	<p>In principle, the brightness decreases with age. This is not a fault: it is normal behaviour.</p> <p>If the brightness decreases from a very early stage, check the following:</p> <ul style="list-style-type: none"> • Automatic dimming due to a temperature problem (flashing display). Allow the light source to cool and ensure that it has cooled adequately before switching back on. • Light guide is faulty and passes less light than previously → replace the light guide • Light guide end-piece or light opening of the light source is dirty → clean them carefully (ensure that light source is cool and not energised!)

Maintenance

No service and maintenance work may be carried out when the device is in use.

This product is essentially maintenance-free. Certain components (e.g. electronics, fans) have a limited service life, resulting in an increase in the probability of a fault developing after a few years' operation. In such cases, the appliance must be repaired at an authorised service centre.

Regularly remove dirt sucked in from the device openings.

At least every 2 years, tests according to EN62353 in the current version are to be carried out and documented. These tests include, as a minimum, protective conductor tests, high-voltage resistance tests, insulation resistance, leakage current and function tests, but are not limited to these. In addition, all tests prescribed by the operator and by law are to be carried out. Check that all electrical and mechanical connections are secure. This work may be performed only by qualified staff. Staff testing and assessing the safety, must be able to recognise possible effects and risks caused by not complying with the specifications of the relevant devices.

Cleaning

Unplug the product from the mains before cleaning and leave it to cool for a sufficient length of time (at least 15 minutes).

Use a dry cloth for cleaning surfaces, or a damp cloth with water and 5% ethyl or isopropyl alcohol; however, under no circumstances should you use combustible or flammable liquids. Commercially-available cleaning and disinfecting wipes may also be used. However, ensure you test on an inconspicuous area in advance, to check that the surface to be cleaned and watermarks do not cause damage. Overly-intensive cleaning or the use of high concentrations can cause scratches, stains or other damage to the surfaces. The light guide must also be cleaned, in particular the area where the light enters and exits.

Remove the light-guide adapter and check the cleanliness of the light-emitting surface (glass surface of 5mm diameter and surrounding parts). If necessary, clean this area gently with a cotton swab, paying attention that no lint is left behind, as these will burn during operation and may damage the device. The device must always be disconnected from the mains before the above process, otherwise there is a risk of serious eye injuries.

If dust has accumulated around the vents, it should be carefully removed. If necessary, use a suitable vacuum cleaner, though you should ensure that the cleaner does not touch the device.

Allow it to dry thoroughly before starting up the device. If you accidentally spill something into the interior of the unit, send the device to the service centre.

Technical data

	F5000 MED
Size (wxdxh)	approx. 295x250x100 excluding projecting parts
Weight	approx. 5kg
Bulb	High Power LED
Brightness	approx. 1000lm with fused fibres of Ø 5x1000mm
Colour temperature	approx. 6900K
Maximum fibre diameter	5mm
Light guide connection	Adaptersystem
Service life L ₅₀	approx. 25,000h (50% of initial brightness)
Operating life	5 years if used as follows: 8h per day 3000h per year 5 years
Cooling	Axial fan
Power supply	100-240VDC 50-60Hz
Power consumption	approx. 90 Watt, approx. 1 Watt during standby
Operating conditions	+10° to +35°C, 20-80% rel. air humidity, 600-1060hPa
Storage and transport	-10° to +50°C, 20-80% rel. air humidity, 600-1060hPa
Rear connections	USB and 2.5mm mono jack for service purposes (covered) Potential equalisation connector
Display	Graphic display for indication of operation
Housing protection class	IP20
Standards	EN60601-1, EN60601-1-2, EN60601-2-18
Medical class	Class 1
EN62471:2008 classification	Risk Group 2
Application part	The application part is the applicator such as an endoscope or microscope, mounted onto the light guide. These applicators are not part of this product. Light guides do not constitute application parts, but the patient may come into contact with them depending on the application (structural insulation in the device of Type CF).
Fuses	2x T2,5AH (250V)

Disposal

Do not dispose of the product in household waste, take it to the appropriate collection and disposal points. Comply with the national and local disposal guidelines and laws. Please help to keep our environment clean.

The development of this product is subject to continuous development. Although we try to keep this operating manual up to date, we reserve the right to change device data at any time without prior notice. The operating manual is produced subject to errors and printing errors.

Edition: 28.03.2018
Doc. No: 633-18-011(b)

Photonic Optische Geräte GmbH&CoKG
Seeböckgasse 59
A-1160 Vienna
 +43-1- 486 56 91-0
 +43-1- 486 56 91-47
 office@photonic.at
www.photonic.at