**▼PHOTONIC**

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<sup>7</sup> Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE Band 7

<sup>8</sup> WLAN 802.11 a/n

Operating Instructions for  
Medical Professionals  
**MINOR SURGICAL LIGHTS**

**ATO ML 600B**  
**ATO ML 1000B**  
**ATO ML 1000K**

(Figure: ATO ML 600B / 1000B)

**MEDICAL PRODUCT**  
Class I according to EU 2017/745 (MDR)

## Dear Customer,

on behalf of PHOTONIC, we hope you enjoy your new ATO ML 600B, ML 1000B or ML 1000K minor surgical lights!

**Attention:** Please read these operating instructions carefully before using the unit for the first time in order to take full advantage of the unit and to avoid damage or injury.

Please use the following order numbers:

<b>Approved accessories (designation)</b>	<b>REF NO.</b>
ATO ML 600 B lamp head	4500-40XX
ATO ML 1000 B lamp head	4500-50XX
ATO ML 1000 K lamp head	4500-30XX
ATO sterilizable handle (part 1/2)	4500.04-020
ATO sterilizable handle (part 2/2)	4500.04-028
ATO ACSwing stand base	4510.0702.XX
ATO ACSwing mobile stand	4510.0701.XX
ATO ACSwing wall mount	4520.0701.XX
ATO ACSwing ceiling mount	4530.0700.XX
ATO ACSwing 1000 ceiling mount	4530.0710.XX
ATO ACSwing 800 ceiling mount	4530.0708.XX
ATO ACSwing 600 ceiling mount	4530.0706.XX
ATO ACSwing 400 ceiling mount	4530.0704.XX
ATO ACSwing 200 ceiling mount	4530.0702.XX

## Table of Contents

<b>Symbols</b> .....	25
<b>Safety Instructions</b> .....	26
<b>Brief Description</b> .....	28
<b>Operation – Positioning</b> .....	29
<b>Operation – Brightness / Endo</b> .....	30
<b>Operation - Color Temperature (ML 1000K)</b> .....	31
<b>Safety Functions / Error Messages</b> .....	32
<b>Preparation of Sterile Accessories</b> .....	33
<b>Disinfection and Cleaning</b> .....	35
<b>Installation and Maintenance</b> .....	36
<b>Disposal</b> .....	37
<b>Incidents and Messages</b> .....	37
<b>Technical Data</b> .....	38
<b>Photometric Data (EN 60601-2-41)</b> .....	39
<b>Electromagnetic Compatibility (EMC)</b> .....	40
<b>Emitted Interference EMC</b> .....	40
<b>EMC Interference Immunity (general)</b> .....	41
<b>EMV Interference Immunity (mobile radio)</b> .....	42

## Symbols

Medical device

Article or serial number

Reference to the operating instructions

**IP42**

Protection against solid foreign bodies and falling dripping water at 15° inclination

Manufacturer / Date of manufacture

Temperature range for transport and storage (from delivery: 15 weeks)

Transport / Storage air pressure range

Transport / Storage humidity

EU - Mark of Conformity

Do not dispose of the product in general domestic waste!

**Please be sure to observe the symbols and notes in the enclosed operating instructions for the supplied support arm system!**

## Safety Instructions

Danger	Note
 Explosion	Do not use in potentially explosive areas, in oxygen-enriched areas or near flammable gases!
 Attraction	The light must not be placed in the vicinity of strong magnetic fields (e.g. MRI system)! Strong attraction!
 Heat	Do not cover the heat sink on the top of the light! Risk of overheating and burns!
 Light	Ensure a maximum total irradiance $< 1000 \text{ W/m}^2$ when using several lights simultaneously!
 Electricity	Only use the enclosed power supply unit and an equipotential bonding conductor! Connect the unit to the mains only after having grounded it with a protective earth conductor!
 Germs	Avoid blocking the displacement ventilation in the operating room! Position the light at an angle to the air flow!
 Failure	The light does not have a fail-safe function in case of power failure! If necessary, use an uninterruptible power supply or have a second unit ready!
<b>The manufacturer is not liable for improper use or repairs of the light!</b>	

Danger	Note
 <b>EMC</b>	<p>Ensure sufficient distance from mobile radio devices. Do not stack any units on top of each other!</p>
 <b>Load</b>	<p>Only use the supplied support arm system and do not attach any additional equipment! Risk of sudden falling of the lamp head!</p>
 <b>Electricity</b>	<p>The mains cable must be unplugged in order to completely disconnect the light from the mains. When using a ceiling mount, it must be possible to switch off the power supply via an on-site switch!</p>
 <b>Tipping</b>	<p>When using the mobile stand, always ensure that the floor is level to avoid accidental tipping over!</p>
 <b>Transport</b>	<p>Only move the mobile light when it is folded up! Disconnect all cables before moving the unit!</p>
 <b>Spring Force</b>	<p>During disassembly of the lamp head, the spring arm can spring open abruptly! Check the connections regularly!</p>
 <b>Maintenance</b>	<p>Assembly and maintenance work may only be performed by trained and PHOTONIC-certified personnel!</p>
 <b>Accessories</b>	<p>For the light, only use accessories that are listed as approved genuine components (see Accessories List / page 2)!</p>

## Brief Description

<b>Intended use</b>
<p>The ML 600B / 1000B / 1000K minor surgical lights are designed for non-invasive, superficial illumination of the entire human body in the visible spectral range for the purpose of examinations, outpatient and inpatient treatments, and surgical procedures by medical professionals. The illumination is used only for optimum visibility of the examination or surgical area and has no diagnostic or therapeutic effect. Illumination is external to the body and the equipment does not come into contact with patients. The ML 600B, ML 1000B and ML 1000K minor surgical lights are intended for use in all rooms used for medical purposes (group 0, 1, 2) – especially operating rooms.</p>
<b>Key Performance Features</b>
<p>Provision of illumination for the surgical areas</p>
<b>Indication</b>
<p>The ML 600B / 1000B / 1000K minor surgical lights are used to illuminate examination or surgical areas to support in diagnoses and treatments, especially during surgical operations.</p>
<b>Contraindication</b>
<p>The lights are not designed for use in dentistry or MRI examinations.</p>
<b>Residual Risks</b>
<p>The light does not have a fail-safe function in case of power failure. In the event of strong collisions, the light may fail and fragments may fall into the operating area.</p>

## Operation – Positioning

The light is positioned primarily by means of the blue, sterilizable handle (G). Non-sterile operators should use the side handle (S) on the light head.

**Note:** During surgery, the blue handles (G) often become non-sterile. Therefore, keep additional handles ready for replacement! (for supplementation: see Accessories list / page 2)

## Operation – Brightness / Endo

Pressing the (|) button (1) turns the light ON at the lowest brightness setting or OFF. The quick start only turns on the light.

**Quick start:** slight rotation of the handle (G)

The brightness of the light can be adjusted by rotating the blue handle (G) or by using the +/- buttons (3,4) on the top of the lamp head.

The set brightness is indicated by the blue indication lights (L) on the handle:


L	max. 60 klx	max. 100 klx
I	10 / 20 klx (ENDO)	10 / 20 klx (ENDO)
II	30 klx	30 / 40 klx
III	40 klx	50 / 60 klx
IIII	50 klx	70 / 80 klx
IIIII	60 klx	90 / 100 klx

## Operation - Color Temperature (ML 1000K)

The light always starts in FULL mode. Pressing the M button (2) switches the light to MIX mode with variable color temperature.

<b>FULL</b> 	4300 K (fixed, 100 klx)	<b>MIX</b> 	3500 - 5000 K (variable, 60 klx)
--	-------------------------------	---	-------------------------------------

The current operating status is indicated via the indication bar on the side of the housing:

Press the M button (2) multiple times to obtain the desired color temperature<sup>1</sup>:

3500 K (MIX I)	3800 K (MIX II)	4100 K (MIX III)
4400 K (MIX IV)	4700 K (MIX V)	5000 K (MIX VI)

<sup>1</sup> Use MIX I or MIX II for optimum red tones ( $R_9 \geq 90$ )

## Safety Functions / Error Messages

The light has built-in safety features to prevent damage to the user and the light, as well as latency after power failures.

### Protection Against Overheating

In the event of overheating, the light switches to an emergency lighting mode with a maximum of 40 klx. The two uppermost indicator lamps on the handle flash. In the event of temperatures exceeding 80 °C at the heat sink, the light automatically switches off the supply voltage and can only be switched on again when the temperature falls below 50 °C.

### Undervoltage

In the event of undervoltage, the light will turn off. All indication lights on the handle flash. If it was on before, the light will return to its previous state once the failure is over.

### Power Outage

In the event of power failure, the light will turn off. If it was on, the light will return to the previous state (brightness, mode) once the failure is over.

### Self-test

The light performs a self-test in the background at regular intervals. If the self-test detects a serious electrical defect, the light is set to a permanent fault state for the protection of the user. All indication lights on the handle flash 5 times when the stand-by button is pressed. In this case, please send the light in for repair!


## Preparation of Sterile Accessories

The light is equipped with a blue, removable and sterilizable handle. This can be disinfected and sterilized by steam up to 134 °C.

To remove the handle (G), press the release buttons (A) on the top and pull down the handle (G) with the release buttons (A):

**Attention:** Before each use (incl. first use) in a sterile environment, the blue handle must be disinfected and sterilized! There is a risk of cross-contamination!

The following procedures have been validated by the manufacturer for reprocessing:

<b>Mechanical Disinfection</b>	
Pre-cleaning	3 min. pre-rinse
Cleaning	5 min at 48 °C, then 2 min cold rinse and rinse with deionized water.
Disinfection	5 min. at 93 °C
Drying	15 min. at 95 °C
Cleaner	Mediclean forte, Dr. Weigert company
Unit	WD 290, Belimed company (program 1)
<b>Sterilization</b>	
Packaging	Blue handle in paper/laminate
Unit	Sterilizer class B, LISA 522 (S/N 08-0794), W&H company
Procedure	Fractionated pre-vacuum process
Temperature	134 °C
Holding time	18 min.

**Note:** In case of deviating processes, it is the responsibility of the user to validate the suitability and effectiveness of the reprocessing!

The handle is designed for up to 1000 reprocessing cycles and must be replaced every 2 years (for supplementation: see accessories list / page 2).

**Attention:** Do not use handles that are damaged or outdated!

## Disinfection and Cleaning

To clean, wipe the surfaces with a slightly damp cloth. If necessary, a mild soap solution can be used. Then wipe dry with a soft anti-static cloth.

**Attention:** Before disinfection and cleaning, the light must be disconnected from the mains!

The light incl. support arm system should be decontaminated on a daily basis or immediately after contamination with a suitable disinfectant (e.g. Meliseptol<sup>1</sup>, Braun) by wipe disinfection.

**Attention:** Wipe cleaning and disinfection only!  
Never spray cleaning agents or disinfectants directly onto the surfaces!

Hygiene guidelines as well as corresponding safety measures for the cleaning and disinfection procedures to be used must be defined by the operator. Please refer to the guidelines of your health care facility!

**Attention:** To avoid the embrittlement of surfaces, do not use agents > 60% alcohol!

---

<sup>1</sup> Recommended and validated by the manufacturer

## Installation and Maintenance

Installation may only be performed by qualified, technical personnel in accordance with the service instructions available from the manufacturer.

**Note:** Your supplier is informed and trained about the scope and content of the installation and maintenance work.

The light incl. support arm system must be inspected annually by service personnel in accordance with EN 62353 and, if necessary, serviced. The inspections performed must be documented and filed by the operator.

**Attention:** Before each use, check the system incl. the handle for visible cracks or deformations as well as for a firm and smooth-running connection between the lamp head and the support arm system!

**Attention:** Loose connections may cause the lamp head to crash unexpectedly!

Do not operate the light if you have any doubts about its mechanical stability, function or electrical safety!

## Disposal

At the end of the product life (10 years), the unit must be disposed of properly and the components separated accordingly.

**Attention:** Do not dispose of the product in general domestic waste!

For proper disposal, contact an approved waste removal company. You can request this from your environmental officer or local government.

**Attention:** Do not disassemble the spring arms and articulated joints! Risk of injury due to pretensioned compression springs!

**Attention:** To avoid environmental contamination, the product must be properly cleaned and disinfected or sterilized before disposal!

## Incidents and Messages

In the event of any serious incident in which the light could be causally involved, the manufacturer and the national health authority must be notified immediately!

## Technical Data

Lamp Head	ML 600B	ML 1000B	ML 1000K
Rated Voltage	24 V DC	24 V DC	24 V DC
Rated Current	1.0 A	1.4 A	1.4 A
Power (elec.)	24 W	33 W	33 W
Protection Class	IP 42	IP 42	IP 42
Weight	2.9 kg	2.9 kg	2.9 kg
Complete system (incl. support arm system and power supply unit)			
Protection Class	I / II*		
Continuous operation	YES (10 years)		
Mains voltage	100 - 240 V AC		
Frequency	50/60 Hz		
Max. power	60 W (power supply)		
Fuses	Ceiling mount	Mobile stand	
Fuse type (IEC 60127)	Primary 250 V T 800 mA H L 5 x 20 mm	Primary 250 V T 800 mA H L 5 x 20 mm	
(* wall mount only: protective conductor is designed as functional grounding).		Secondary 250 V M 2000 mA L 5 x 20 mm	
Ambient Conditions - Operation			
Ambient temperature		10 – 40 °C	
Rel. humidity (non-condensing)		30 – 75 %	
Air pressure		700 – 1060 hPa	
Ambient Conditions - Transport <sup>1</sup>			
Ambient temperature		-25 – 70 °C	
Rel. humidity (non-condensing)		10 – 75 %	
Air pressure		500 – 1060 hPa	

<sup>1</sup> Storage conditions up to 15 weeks from delivery date

## Photometric Data (EN 60601-2-41)

<b>ML 600B<sup>A</sup>   1000B<sup>B</sup>   1000K<sup>C</sup> (in 100 cm, mean values)</b>	
Central illuminance [klx]	60 <sup>A</sup> /100 <sup>B,C</sup>
Total irradiance [W/m <sup>2</sup> ]	220 <sup>A</sup> /353 <sup>B</sup> /390 <sup>C</sup>
Illuminated field diameter D <sub>50</sub> [mm]	140
Illuminated field diameter D <sub>10</sub> [mm]	210
D <sub>50</sub> / D <sub>10</sub> – Distribution [-]	0.66
Illumination depth 60% [mm]	≥ 1100
Minimum working distance [cm]	70
Color rendering index R <sub>a</sub> [± 10%]	≥ 95
Color rendering index R <sub>9</sub> [± 10%]	≥ 93 <sup>A</sup> /84 <sup>B</sup> /79 <sup>C</sup>
Shading without screen [%]	100
Shading with one screen [%]	0
Shading with two screens [%]	66 <sup>A</sup> /63 <sup>B</sup> /62 <sup>C</sup>
Shading with tube [%]	100
Shading with tube/screen [%]	66 <sup>A</sup> /63 <sup>B</sup> /62 <sup>C</sup>
LED lifetime [h]	≥ 50,000
<b>ATO ML 600B - Color temperature / brightness</b>	
Color temperature [K]	4000
Brightness adjustment [klx]	10 – 60
<b>ATO ML 1000B - Color temperature / brightness</b>	
Color temperature [K]	4000
Brightness adjustment [klx]	10 – 100
<b>ATO ML 1000K - Color temperature / brightness</b>	
Color temperature [K, ± 200] (variable, MIX - mode)	3500 – 5000 (in 6 steps)
Brightness adjustment (MIX) [klx]	10 - 60
Color temperature [K, ± 200] (FULL - mode)	4300
Brightness adjustment [FULL]	10 - 100

## Electromagnetic Compatibility (EMC)

Electromagnetic fields may interfere with the light. Likewise, the light may affect other medical units and impair their function. Therefore, observe the limit values specified here!

**Attention:** These limit values only apply to the lights in combination with the approved support arm systems and accessories (see page 2). The use of other or additional components is strictly prohibited!

## Emitted Interference EMC

Please note that the light may only be operated in an electromagnetic environment as specified below (see table)!

Phenomenon	Standard section	Test level <sup>1</sup>
Emitted interference <sup>2</sup>	CISPR 15	Group 1, Class B
Distortion <sup>3</sup>	IEC 61000-3-2	Class A
Flicker <sup>4</sup>	IEC 61000-3-3	compliant

---

<sup>1</sup> Environment in home health care settings (including professional health care facilities).

<sup>2</sup> Power guided and radiated interference emissions

<sup>3</sup> Harmonic distortion

<sup>4</sup> Voltage fluctuations and flicker

## EMC Interference Immunity (general)

Please note that the light may only be operated in an electromagnetic environment as specified below (see table)!

Phenomenon	Standard section	Test level <sup>1</sup>
Discharge <sup>2</sup>	IEC 61000-4-2	Contact: $\pm 8$ kV Air: $\pm 2, 4, 8, 15$ kV
Burst <sup>3</sup>	IEC 61000-4-5	$\pm 2$ kV @ 100 kHz
Surge cable <sup>4</sup>	IEC 61000-4-5	$\pm 1$ kV
Surge ground <sup>5</sup>	IEC 61000-4-5	$\pm 2$ kV
Dip <sup>6</sup>	IEC 61000-4-11	0% UT: 1 period 70% UT: 25/30 per.
Magnetic field <sup>7</sup>	IEC 61000-4-8	30 A/m
HF fields I <sup>8</sup>	IEC 61000-4-3	10 V/m @ 0.08 – 2.7 GHz
HF fields II <sup>9</sup>	IEC 61000-4-6	10 V <sub>RMS</sub> @ 0.15 – 80 MHz

If you have questions, contact your service technician!

---

<sup>1</sup> Environment in home health care settings (including professional health care facilities).

<sup>2</sup> Static electricity discharge

<sup>3</sup> Fast transient electrical disturbances

<sup>4</sup> Surge voltage cable to cable

<sup>5</sup> Surge voltage cable to ground

<sup>6</sup> Voltage dips

<sup>7</sup> Magnetic fields with energy technology rated frequencies

<sup>8</sup> High-frequency electromagnetic fields

<sup>9</sup> Conducted interference induced by high-frequency fields

## EMV Interference Immunity (mobile radio)

Please ensure sufficient distance between the light and mobile devices to avoid electromagnetic interference!

Observe the following limit values (see table)!

Frequency [MHz]	Modulation	$P_{\max}$ [W]	d [m]	$E_p$ [V/m]
380 - 390 <sup>1</sup>	PM 18 Hz	1.8	0.3	27
430 - 470 <sup>2</sup>	FM 5/1 kHz <sup>3</sup>	2.0	0.3	28
704 - 787 <sup>4</sup>	PM 217 Hz	0.2	0.3	9
800 - 960 <sup>5</sup>	PM 18 Hz	2.0	0.3	28
1700 - 1990 <sup>6</sup>	PM 217 Hz	2.0	0.3	28
2400 - 2570 <sup>7</sup>	PM 217 Hz	2.0	0.3	28
5100 - 5800 <sup>8</sup>	PM 217 Hz	0.2	0.3	9

### Key:

PM	Pulse modulation
FM	Frequency modulation
$P_{\max}$	maximum power of the radio
d	Distance to the radio
$E_p$	Immunity test level

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Photonic Optische Geräte GmbH & Co KG