

FIBRE-OPTIC LED LIGHT SOURCE F1/F1R

COLD LIGHT SOURCE	LED F1 (633-20-00299)	LED F1 R (633-20-00399)
Lamp	High-Power LED	
Control range	min. 30:1	
Colour temperature	5 500 K (daylight)	
Power input	22 W	
Power supply	24 V DC +/- 5 %, 1000 mA, socket for hollow pin 5.5 x 2.1 mm	
Brightness control	infinitely adjustable using rotary potentiometer; with index marks	
LED service life	50 000 h	
Light guide connection	Photonic standard (Ø 15 mm, galvanic isolation)	
Remote socket	no	yes
ESD socket	for 4 mm banana plug for potential equalisation	
Power cord	plug-in	
Dimensions (L x W x H)	96 x 110 x 175 mm without controls	
Weight	2 500 g	
Type of protection (protection class)	III	
Conformity	CE (EMC Directive 89/336/EEC), RoHS	